




MICROPLANET LABORATORIOS, SL

www.microplanet-psl.com

info@microplanet-psl.com

Tel. 93 590 62 92

PRODUCTOS Y SERVICIOS PARA LABORATORIOS





Seguridad alimentaria




Mejorando la calidad y la productividad en su trabajo

Systems for microbial identification directly from food and surfaces

Description	Packaging	Ref.
Food System System for detection and identification of pathogenic microorganisms from food. VALIDATED to ISO 16140 standard Microbiology of food and animal feeding stuffs - Protocol for the validation of alternative methods	40 test 8 test	71680 79680
		
HACCP System Plus System for enumeration and identification of pathogenic microorganisms from surfaces.	40 test 8 test	71655 79655
		

Coliforms in water testing

Description	Packaging	Ref.
		
Water Detect Coliforms Chromogenic medium for detecting <i>E. coli</i> and others coliforms in water samples.	20 tubes 12 bottles x 100 mL	630600 630601

Pesticide residue detection

Description	Packaging	Ref.
		
MAIA Pesticide MultiTest Colorimetric system for the detection of organophosphate, organochloride and carbamate pesticides residues in water, food and drink.	96 test	79700

Antibiotic residue detection		
Description	Packaging	Ref.
MeRA Test Detection of antibiotic and sulphamide residues in meat by <i>Geobacillus stearothermophilus</i> .	50 vials	80356
MiRA Test Detection of antibiotic and sulphamide residues in milk by <i>Geobacillus stearothermophilus</i> .	50 vials	80355
PAR TEST Agar Medium with spores of <i>Geobacillus stearothermophilus</i> subsp. <i>calidolactis</i> for residual antibiotic detection in biological liquids (milk, blood and urine).	20 plates 100 plates	11200 11200*
BLANK DISCS 6 mm Paper discs.	250 discs	9999
Antibiotic Test Antibiotic residue detection in urine and other biological fluids by <i>Bacillus subtilis</i> .	50 discs	80350
Rapid Antibiotic Test Antibiotic residue detection in foodstuffs by <i>Geobacillus stearothermophilus</i> var. <i>calidolactis</i> .	50 tests	80351
Suspension Spore BAS E8 Sterile Suspension of <i>Bacillus subtilis</i> spores in distilled water for antibiotic residue detection in milk and other biological fluids (1-5x10 ⁸ CFU/mL).	10 mL 2 mL	82012 82212
Suspension Spore BAS E7 Sterile Suspension of <i>Bacillus subtilis</i> spores in distilled water for antibiotic residue detection in milk and other biological fluids (1-5x10 ⁷ CFU/mL).	10 mL	82013
Suspension Spore BAS E6 Sterile Suspension of <i>Bacillus subtilis</i> spores in distilled water for antibiotic residue detection in milk and other biological fluids (1-5x10 ⁶ CFU/mL).	10 mL	82014
Suspension Spore BAS E5 Sterile suspension of <i>Bacillus subtilis</i> spores in distilled water for antibiotic residue detection in milk and other biological fluids (1-5x10 ⁵ CFU/mL).	10 mL	82015


Solutions for sample collection		
Description	Packaging	Ref.
SWAB Solution SKELETON Solution to hydrate samples collection swabs (2001/471/CE).	20 vials	20157
SWAB Solution SURFACES Solution to hydrate samples collection swabs before surfaces checking.	20 vials	20160

Swabs w/ and w/o transport medium		
Description	Packaging	Ref.
BIO – CULT AMIES + charcoal Swabs for collection and transport of clinical, environmental samples and microorganisms.	150 Swabs	85307
BIO – CULT AMIES Swabs for collection and transport of clinical, environmental samples and microorganisms.	100 Swabs	85317
BIO CULT CARY BLAIR Swabs for collection and transport of clinical, environmental samples and microorganisms.	150 Swabs	85310
BIO – CULT Stuart + charcoal Swabs for collection and transport of clinical, environmental samples and microorganisms.	150 Swabs	85306
BIO – CULT Stuart Swabs for collection and transport of clinical, environmental samples and microorganisms.	100 Swabs	85308
Sterile Swab in Envelope Dry swabs for clinical and environmental samples collection.	100 Swabs	85270
Sterile Swab without Envelope Swabs for clinical and environmental samples collection.	100 Swabs	85275


Transport medium for Viruses, Chlamydia, Mycoplasma and Ureaplasma		
Description	Packaging	Ref.
VTM	100 tubes x 3 mL	26490
VTM	6 bottles x 200 mL	412450






Organic residue detection on surfaces

Description	Packaging	Ref.
 <p>Clean Test Swab for detecting protein, fat and sugar residues for the hygiene control of surfaces.</p>	50 test	86001
	3 test	86001/3

Surface swabs for microbial detection

Description	Packaging	Ref.
 <p>Contam Swab Listeria</p> <p>Contam Swab E.coli-coliforms</p> <p>Contam Swab Salmonella</p> <p>Contam Swab Contamination</p> <p>Contam Swab MRSA</p> <p>Contam Swab Yeasts & Moulds</p>	30 test	86102
	30 test	86101
	30 test	86103
	30 test	86100
	30 test	86104 #
	30 test	86105 #

ESC (Easy Surface Checking) swab (ISO 18593)

Description	Packaging	Ref.
   <p>ESC Swab Neutralizing Rinse Solution 10 mL</p> <p>ESC Swab Neutralizing Rinse Solution 5 mL</p> <p>ESC Swab Buffered Peptone Water 10 mL</p> <p>ESC Swab Buffered Peptone Water 5 mL</p> <p>ESC Swab D/E Neutralizing Broth 10 mL</p> <p>ESC Swab D/E Neutralizing Broth 5 mL</p> <p>ESC Swab Letheen Broth 10 mL</p> <p>ESC Swab Letheen Broth 5 mL</p> <p>ESC Swab Letheen Broth 1 mL</p> <p>ESC Swab Maximum Recovery Diluent 10 mL</p> <p>ESC Swab Maximum Recovery Diluent 5 mL</p> <p>ESC swab Mycoplasma Transport Broth</p> <p>ESC swab Demi Fraser Broth 10 mL</p> <p>Sampling Template 10 cm x 10 cm sterile</p>	100 test	85601
	100 test	85602
	100 test	85603
	100 test	85604
	100 test	85605
	100 test	85606
	100 test	85607
	100 test	85608
	100 test	85621
	100 test	85609
	100 test	85610
	100 test	85612
	100 test	85614
	60 units	96762

Liofilchem® s.r.l.

Via Scozia, Zona Industriale
64026, Roseto degli Abruzzi (Te)
Italy
tel +39 0858930745
fax +39 0858930330
www.liofilchem.com



MICROPLANET LABORATORIOS, SL
www.microplanet-psl.com
info@microplanet-psl.com
Tel. 93 590 62 92