

## ESC Swab

Torunda en tubo, con tapón de rosca y una amplia gama de diluyentes



**ESC (Easy Surface Checking) Swab**, es una nueva gama de hisopos para muestreo de superficies, especialmente recomendado para la industria alimentaria y farmacéutica.



**Características del producto:**

- Torunda de rayon
- Diluyente incorporado
- Amplia gama de diluyentes disponibles
- Tapón a rosca y hermético a prueba de fugas
- Transporte y almacenamiento a temperatura ambiente
- En cumplimiento de las normas ISO 18593 e ISO 17604



Ref.	Descripción	Envase
2360085601	Neutralizing Rinse Solution	100 x 10 mL
2360085602	Neutralizing Rinse Solution	100 x 5 mL
2360085603	Buffered Peptone Water	100 x 10 mL
2360085604	Buffered Peptone Water	100 x 5 mL
2360085605	D/E Neutralizing Broth	100 x 10 mL
2360085606	D/E Neutralizing Broth	100 x 5 mL
2360085607	Lethen Broth	100 x 10 mL
2360085608	Lethen Broth	100 x 5 mL
2360085609	Maximum Recovery Diluent	100 x 10 mL
2360085610	Maximum Recovery Diluent	100 x 5 mL

**Plantilla de Muestreo 10 x 10**

**Plantilla de Muestreo 10 x 10:**

Dispositivo de plástico estéril utilizado para delimitar una superficie de 100 cm<sup>2</sup>, según procedimiento ISO 18593:2004

Ref.	Descripción	Envase
2360096762	Plantilla de muestreo	60 unid.



**Referencias:**

1. ISO 18593:2004 *Microbiology of food and animal feeding stuffs - Horizontal methods for sampling techniques from surfaces using contact plates and swabs.* International Standards Organisation, Geneva
2. ISO 17604:2003 *Microbiology of food and animal feeding stuffs - Carcass sampling for microbiological analysis.* International Standards Organisation, Geneva
3. ISO 6887-1:1999 *Microbiology of food and animal feeding stuffs - Preparation of test samples, initial suspension and decimal dilution for microbiological examination.* International Standards Organisation, Geneva

